

# VETERINARY MICROSCOPE

## SET UP CHECKLIST

Setting up a microscope involves a series of steps to ensure it's ready for accurate observations. Here's a comprehensive guide:

### 1. Clean the Workspace

- Begin by clearing the workspace of any clutter or debris. Ensure there's enough room for the microscope and necessary accessories.

### 2. Inspect the Microscope

- Visually inspect the microscope for any signs of damage or issues from the previous day. Check for loose knobs, cords, or components.

### 3. Check the Electrical Connections

- Ensure that the microscope is securely plugged into a power source. Verify that the power cord is in good condition without fraying or exposed wires.

### 4. Turn On the Microscope

- Switch on the microscope using the power button. Allow it a few moments to initialize.

### 5. Adjust the Illumination

- Set the illumination to an appropriate level for the type of observations you'll be conducting. Avoid excessive brightness, which can cause glare.

### 6. Clean the Lenses

- Use lens paper and a gentle optical cleaning solution to clean the objective and eyepiece lenses. Start from the center and move outward in a circular motion.

### 7. Check the Light Source

- Inspect the light source to ensure it's working properly. Replace any bulbs that have burned out.



# VETERINARY MICROSCOPE

## SET UP CHECKLIST

### 8. Position the Microscope

- Place the microscope on a stable, level surface. Ensure it's positioned in a way that allows comfortable viewing and easy access to the stage controls.

### 9. Power Up Additional Components

- If your microscope is connected to a camera or computer for imaging, ensure these components are powered up and functioning.

### 13. Perform a Functionality Check

- Test all the microscope's functions, including focusing, stage movement, and any special features specific to your model.

### 14. Document Any Issues

- If you encounter any problems during the set-up process, document them for further evaluation or maintenance.

**By following this set-up routine, you'll ensure that your microscope is in optimal condition for accurate and reliable observations throughout the day and avoids the risk of an issue being found when the microscope is needed for use.**

## GENERAL TIPS

- Always handle the microscope with clean, dry hands to prevent contamination.
- Store the microscope in a clean, dust-free environment.
- Follow the manufacturer's guidelines for specific maintenance recommendations and avoid using unauthorised accessories or lubricants.
- By following this maintenance guide, you can extend the lifespan of your veterinary microscope and ensure it consistently delivers accurate and reliable results in your diagnostic procedures.



# MICROSCOPE DO'S AND DON'TS



## DO

### Handle with Care

- Handle the microscope gently and with clean, dry hands to prevent damage or contamination.

### Clean Lenses Properly

- Use lens paper and a gentle optical cleaning solution to clean the objective and eyepiece lenses. Start from the center and move outward in a circular motion.

### Use Lens Caps

- Always replace lens caps when the microscope is not in use to protect the lenses from dust and scratches.

### Start with Low Magnification

- Begin your observations with the lowest magnification objective to locate and center your sample. Gradually increase magnification as needed.

### Adjust Illumination

- Use the microscope's light source adjuster to control brightness. Avoid excessive light, as it can cause glare and reduce image quality.

### Use Coarse and Fine Focus Knobs Properly

- Start with the coarse focus knob to bring the specimen into rough focus, then use the fine focus knob for precise adjustments.

### Keep Slides Clean

- Ensure slides are free from debris, fingerprints, or stains before placing them on the microscope stage.

### Store Properly

- Cover the microscope with a dust cover when not in use. Store it in a clean, dry environment away from potential sources of damage.

### Follow Manufacturer's Instructions

- Adhere to the specific guidelines and recommendations provided by the microscope's manufacturer.

## DON'T



### Force Knobs or Adjustments

- Avoid applying excessive force when using focus knobs or making adjustments. This can damage the microscope's internal mechanisms.

### Touch Lenses with Fingers

- Refrain from touching the objective or eyepiece lenses with your fingers, as oils and dirt can degrade image quality.

### Use Abrasive Materials for Cleaning

- Avoid using rough materials or strong chemicals for cleaning lenses, as they can scratch or damage the glass.

### Adjust Lighting Drastically

- Do not suddenly increase the brightness to maximum. Gradual adjustments are more effective and prevent overexposure.

### Skip Proper Calibration

- Ensure that your microscope is properly calibrated for accurate measurements. Neglecting this step can lead to inaccurate data.

### Store in Humid or Dusty Environments

- Avoid storing the microscope in areas with high humidity or excessive dust, as these conditions can cause damage over time.

### Lean on the Microscope

- Do not lean on or apply pressure to the microscope. This can misalign components or potentially cause structural damage.

### Disassemble without Training

- Do not attempt to disassemble or repair the microscope unless you are trained and authorised to do so. Seek professional assistance if needed.

**Following these dos and don'ts will help ensure the longevity and optimal performance of your microscope for accurate and reliable observations.**



# **VETERINARY** MAINTENANCE **MICROSCOPE** AGENDA

## **DAILY MAINTENANCE**

### **1. Cleaning**

- Use a soft, lint-free cloth to wipe down the microscope body, focusing knobs, and objectives to remove dust and fingerprints.
- Use lens paper and a gentle optical cleaning solution to clean the objective lenses. Avoid using harsh chemicals or solvents.

### **2. Checking for Loose Parts**

- Inspect for any loose knobs, screws, or components. Tighten them gently if necessary.

### **3. Inspecting Electrical Components**

- Ensure that the power cord is in good condition with no fraying or exposed wires. Check the plug for any signs of damage.

## **WEEKLY MAINTENANCE**

### **1. Condenser and Light Source**

- Clean the condenser lens and light source with a soft brush or compressed air to remove dust or debris.

### **2. Stage and Mechanical Parts**

- Inspect the stage for smooth movement. Lubricate the mechanical parts if needed with a manufacturer-recommended lubricant.





# **VETERINARY** MAINTENANCE **MICROSCOPE** AGENDA

## **MONTHLY MAINTENANCE**

### **1. Focus Mechanism**

- Check the focus knobs for smooth operation. Lubricate if necessary, following manufacturer guidelines.

### **2. Filters**

- If your microscope has filters, inspect them for dirt or scratches. Clean or replace them as needed.

### **3. Objective Lenses**

- Remove the objectives and clean them using a gentle optical cleaning solution and lens paper. Ensure they are dry before reattaching.

## **ANNUALLY OR AS NEEDED**

### **1. Professional Servicing**

- Arrange for a professional technician to perform a comprehensive maintenance and calibration service. This may include aligning optics, checking electrical connections, and replacing worn components.

### **2. Calibration**

- Ensure that the microscope is calibrated for accurate measurements. This may involve using a calibration slide or having a technician perform the calibration.

### **3. Documentation**

- Keep a log of all maintenance activities, including dates, tasks performed, and any issues or repairs addressed and fixed. This helps track the microscope's history and informs future maintenance needs.

