
Floline V5 Infusion Pump with Built In Fluid Warmer







Flexible Adjustment:

30~40°C continuously adjustable, temperature control accuracy ±1.0°C.



Multiple Protections:

The heating will be automatically turned off when the temperature exceeds 50°C to ensure the safety of infusion.



Rapid Heating:

4 sets of heating units quickly reach the set temperature.

Extended Heating Tube:

1100mm long heating tube, which can effectively heat the medicinal liquid and reduce the loss of drug activity.

Safe and Reliable:

Features a flexible three chamber structure, fully sealed with a waterproof design.

Precise Temperature Control: 2 NTC temperature detection to ensure more accurate heating control.

Fluid Warmer works exclusively with Floline V5 Infusion Pump



Check out our website at www.vet1.com.au or e-mail sales@vet1.com.au or call us on **1300 378 713** to discuss anything Veterinary Equipment related. We are Committed to Quality.